



**Service information**

**Reference no.: SDA\_113343**

**Symptom**

Single components of Flow Selector Faucet not available.



**Cure**

**Content:** The disassembly/assembly of flow selector faucet may cause:

- incorrect calibration
- water leakage
- incorrect functionality of flow selector faucet

**Cure:** To prevent issues due to incorrect assembly the flow selector faucet is just available as assembled part.

**Actions:** Order the Flow Selector Faucet Assembly:

<p><b>For Intelia, Intuita, Gaggia Naviglio, OCS Vienna, Moltio, Philips 2100: 17000114 996530068077 FLOW SELECTOR FAUCET NPR</b></p>	
<p><b>For Xsmall Vapore, Phi2000 Vapore, Gaggia Brera: 11012381 996530006383 FLOW SELECTOR FAUCET V2 P0057.</b></p>	

**Do not use cannibalized parts from other machine.**

In case you have stocks of the singular components:

- 12000893 996530007948 OR ORM 0120-20 BLACK EPDM;
- 11008104 996530004722 MUSHROOM VALVE FOR FLOW SELECTOR P0057;
- 11008105 996530004723 FLOW SELECTOR VALVE SUPPORT P0057;
- 11008077 996530004709 STEEL FLOW SELECTOR VALVE SPRING P0057;
- 11008106 996530004724 FLOW SELECTOR VALVE CAP P0057.

do not use these any more for repair but contact the spareparts order desk of SAEKO in Italy for the appropriate action.

**ADDITIONAL INFORMATION:**

**SERIAL NUMBER** (First serial where solution has been implemented in production):

None.

**IPM NUMBER** (Philips Integrate Problem Management documentation reference number):

None.

**INFORMATION FOR CONSUMER:**

None.

**HISTORY:**

26/04/2016 - First Release